


Aluminum Alloy Folding Stretcher

I. Technical Parameters:

1. Product Size (L×W×H): 206×52×13cm
2. Folding Size (L×W×H): 103×18×9cm
3. Carton Size (2pcs) :: 108×18.5×18cm
4. N.W: 4.7kg G.W.: 11kg
5. Loading Bearing: ≤159kg



II. How to use:

The folding stretcher with outrigger is made of high strength aluminum alloy material and the waterproof leather materials. This stretcher with outrigger can be folded in length and width, easy for transportation and storage. The orange fabric is fire proof, water proof, anti-cracking and easy to clean.

III. Applied for:

Approved by clinic practice, this item is suitable for ambulance, hospital battlefield gym, to carry the wounded & patients. It also can be used in general operation.

IV. Cautions:

1. When carrying patients with this stretcher, be sure lock tightly fasten the safety belt for safety sake.
2. In general operation, be careful not to out torn the stretcher surface.

V. Maintenance:

1. Keep clean as routine (incl. sterilization).
2. Often check if any parts loose or not.

VI. Storage & transportation:

1. Store this product in the place of damp proof & non-corrosion.
2. Common transporting vehicle can carry this product.